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Stump shrinkers



---

Gebrauchsanweisung · Instructions for use · Mode d'emploi · Istruzioni d'uso · Instrucoes de utilizacao · Instrucciones de uso · Gebruiksaanwijzing · Brugsanvisning · Bruksanvisning · Kullanım talimatları · Upute za uporabu · Инструкция по эксплуатации · ارشادات · Navodila za uporabo

---

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## Intended purpose

Compression therapy comprises measures of external, surface pressure application and is used in cases of damage to the venous and/or lymphatic vessel function/structure on the amputation stump with swelling, oedema, possibly changes to the skin, for oedema reduction and prevention of stump volume fluctuations. The aim of surface pressure is to prevent the formation of oedema and to support venous return and lymph drainage. Compression garments exert a perfect pressure gradient that gradually decreases from the bottom to the top. It involves the use of compressive and elastic double-stretch fabrics or knitted fabrics that wrap right around the extremities. As a general rule, compression therapy is a long-term treatment requiring appropriate products that are used by the patients themselves when medically indicated and as prescribed by a doctor.

## Notes for use

Your specialist retailer will be happy to show you the correct way to use the product. Wear your compression garment as prescribed by your doctor. Put on your compression garment while you're sitting down to prevent injuries or falling. Rough, unfiled finger nails or sharp-edged rings etc. may damage the very fine knitted fabric. We therefore recommend the use of **Juzo special gloves**.

1. Turn the shrinker inside out, except for its end part.
2. Take the shrinker, folded over itself, in both hands, and slip your stump into it.
3. Now pull up the shrinker evenly and free of wrinkles – without twisting it or pulling on the top band – and use the flat of your hand to shape it to your leg. If you pull the shrinker up too far, simply slide it back down and slip it up your leg again with less tension this time. Tip: Fold over the silicone border (if present) before putting on the shrinker.

## Please note:

Only wear compression garments if they have been prescribed by a doctor. Compression garments should only be dispensed by properly trained staff at a medical retailer. Your Juzo compression garments are largely resistant to oils, ointments, skin moisture and other environmental influences. The durability of the knitted fabric may be affected by certain combinations of the above. This is why it is advisable to have your compression garments checked regularly at a medical retailer. For stockings with a silicone border, please note that legs should be depilated and free from oils, powder, creams and ointments. Please contact your specialist retailer in the event of damage. Do not repair the product yourself – this may affect its quality and medical effectiveness. For hygiene reasons, this product is intended for one person only.

## Washing and care instructions

Please note the care instructions on the textile label sewn into your compression garment or the printed information. Do up any fasteners (if there are any) before washing! We recommend turning compression garments inside out before washing. Wash compression garments separately when washing for the first time (they may lose colour). If possible, your compression garments should be washed every day on a gentle or delicate wash cycle (40°C). Tip: a laundry net provides additional protection for the knitted fabric. Please do not use fabric softener! When hand washing, rinse well and do not wring. We recommend the gentle Juzo special detergent. You can reduce drying time by placing the compression garment on a thick towel after washing, tightly rolling the towel up and firmly pressing out excess moisture. Then hang up the compression garment loosely. Do not leave the product in the towel and do not dry it on a radiator or in the sun. Please use the delicates setting when drying in a tumble dryer (suitable for tumble-drying: see textile label). Juzo compression garments are not suitable for dry cleaning.

## Material composition

For precise details, please see the textile label sewn into your compression garment.

### Storage information and shelf life

Store in a dry place and protect from sunlight. Standard products have a maximum shelf life of 36 months. The compression garment then has a usage period of 6 months. The usage period depends on correct handling (e.g. care, donning and doffing) and normal use of the products. The usage period of the medical product is printed on the box with an hourglass symbol. Custom-made products are for immediate use. Body measurements may change based on the patient's condition and individual circumstances. Custom-made products are therefore designed for a usage period of 6 months. After the 6-month usage period, body measurements must be re-checked by a medical specialist. In the event of significant changes in measurements of the body part that is being treated (e.g. based on therapeutic success, weight change), re-measurement and supply of a new product may be required earlier, before the 6-month usage period comes to an end.

### Indications

- In cases of damage to the venous and/or lymphatic vessel function/structure on the amputation stump with swelling, oedema, possibly with changes to the skin
- For oedema reduction
- Preventing stump volume fluctuations after removal of prosthetics

Following burns and for scar treatment, we recommend using our special scar therapy products.

### Absolute contraindications

Peripheral artery disease (PAD) stage III/IV, decompensated heart failure, septic phlebitis, phlegmasia cerulea dolens, acute erysipelas, acute bacterial, viral or allergic inflammation, untreated or advanced skin diseases/contact dermatitis, untreated open wounds, swelling of the extremity of unknown cause, untreated malignant tumours.

### Relative contraindications

Pronounced weeping dermatoses, malignant lymphoedema, peripheral artery disease (PAD) stage I/II, peripheral neuropathy in association with diabetes mellitus, sensory disorders of the extremities, compensated heart failure, skin infections, skin diseases (e.g. blistering dermatoses), intolerance to pressure or the ingredients of the product, chronic rheumatoid arthritis. If there are open wounds in the area to be treated, the affected area must be covered with an appropriate wound dressing or bandage before the compression stocking can be applied. No liability can be accepted by the Julius Zorn GmbH Group in the event of failure to comply with these contraindications.

### Side effects

There are no known side effects when used as directed. However, if negative alterations (such as for example skin irritations) should occur while the use of our products has been prescribed, please go to your doctor or your specialist dealer immediately. If an incompatibility should be known against one or several ingredients of this product, please check back with your doctor before use. The manufacturer and any of its affiliates is not liable for damages / injuries which have been caused by improper handling or misuse.

In the case of complaints related to this product, such as damage to the knitted fabric or problems with the fit, please contact the medical retailer directly. Only serious incidents that may lead to a major deterioration of the patient's medical condition or to death should be reported to the manufacturer and the relevant authority in the Member State. Serious incidents are defined in Article 2 (65) Regulation (EU) 2017/745 (MDR).

### Disposal

You can dispose of your compression garments in the usual way. There are no specific disposal instructions.

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**Julius Zorn GmbH** · Juliusplatz 1 · 86551 Aichach · Germany