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





Purpose

The consistent application of uniform surface pressure by the compression garments has a beneficial effect on scar maturation and development. Therapeutic application of compression garments results in decreased blood circulation and increased tissue elasticity. The aim is to reduce the risk of the formation of hypertrophic scars and/or keloids or achieve a reduction in their appearance. The objective of compression therapy is to flatten the scar, increase tissue elasticity and reduce scar rigidity.

Notes for application

Your specialist dealer or therapist will show you how to use it correctly. Wear your compression garments for 23 hours non-stop (or as advised by your doctor). Compression garments are perfect for whole body treatments and they offer a certified UPF 80 UV protection factor in accordance with UV Standard 801. The UV protection only applies to those body regions that are covered with compression garments. Rough, unfilled nails or sharp-edged rings etc. may damage the very fine knitted fabric. Therefore please use Juzo Special Gloves. For compression garments with a silicone border, please note: The relevant part of the body should be depilated and free from oils, powder, creams and ointments.

Compression sleeves with knitted edges or silicone borders

1. Starting from the top, down to half way, turn the sleeve inside out.
2. Pull the sleeve over your hand and work your way up your arm to the elbow.
3. Work the folded section of the sleeve over your elbow and pull it step by step towards your shoulder.
4. Spread out the sleeve's knitted fabric so that there are no creases. When the sleeve is in the right position and fits perfectly, it will neither cause uncomfortable pressure or pinching, nor restrict your movements. Tip: Before donning, fold over the silicone border (if applicable).

Compression sleeves with shoulder cap and strap or bra fastening

5. Place one part of the shoulder strap across your chest and the other over your back. Next, fasten both ends together. Please bear in mind that you can change the length of the shoulder strap using the adjuster if it is too loose or too tight.

Please note:

Only wear compression garments if they have been prescribed by a doctor. Compression garments should only be dispensed

by properly trained staff at a medical retailer. Progressing the success of treatment or a change in condition means that the compression garments will have to be remeasured and, if necessary, new items will need to be made. Your compression garments are largely resistant to oils, ointments, skin moisture and other environmental influences. The durability of the knitted fabric may be affected by certain combinations of the above. This is why it is advisable to have your compression garments checked regularly at a medical retailer. Please contact your specialist retailer in the event of damage. Do not repair the product yourself – this may affect its quality and medical effectiveness. For hygiene reasons, this product is intended for one person only.

Washing and care instructions

Please note the care instructions on the textile label sewn into your compression garment. Do up any fasteners before washing! We recommend turning compression stockings with silicone borders inside out before washing. Wash compression garments separately when washing for the first time (they may lose colour). Your compression garments should be washed every day on a gentle or delicate wash cycle (40 °C). Tip: A laundry net provides additional protection for the knitted fabric. Please do not use fabric softener. When hand washing, rinse well and do not wring. We recommend the gentle Juzo Special Detergent. You can reduce drying time by placing the compression garment on a thick towel after washing, tightly rolling the towel up and firmly pressing out excess moisture. Then hang up the compression garment loosely. Do not leave it in the towel and do not dry on a radiator or in the sun. Compression garments are not suitable for dry cleaning or tumble drying.

Material composition

For precise details, please see the textile label that is sewn into your compression garment.

Storage information and shelf life

Store in a dry place and protect from sunlight. Custom-made products are for immediate use. Body measurements may change based on the patient's condition and individual circumstances. Custom-made products are therefore designed for a usage period of 6 months. After the 6-month usage period, body measurements must be re-checked by a medical specialist. In the event of significant changes in measurements of the body part that is being treated (e.g. based on therapeutic success, weight change), re-measurement and supply of a new product may be required earlier, before the 6-month usage period comes to an end.

Indications

Damage to the skin and healing for skin with hypertrophic scars and keloids. After injuries/burns/operations, to improve the scar/wound healing through to completely healing the scar, to stabilise the surgical outcome.

Absolute contraindications

Peripheral artery disease (PAD) stage III/IV, decompensated heart failure, septic phlebitis, phlegmasia cerulea dolens, acute erysipelas, acute bacterial, viral or allergic inflammation, untreated or advanced skin diseases/contact dermatitis, untreated open wounds, swelling of the extremity of unknown cause, untreated malignant tumours.

Relative contraindications

Pronounced weeping dermatoses, malignant lymphoedema, peripheral artery disease (PAD) stage I/II, peripheral neuropathy in association with diabetes mellitus, sensory disorders of the extremities, compensated heart failure, skin infections, skin diseases (e.g. blistering dermatoses), intolerance to pressure or the ingredients of the product, rheumatoid arthritis. If there are open wounds in the area to be treated, the affected area must be covered with an appropriate wound dressing or bandage before the compression garment can be applied. For chronic lymphoedema or lipoedema, we recommend our lymphological products.

Liability for failure to comply with these contraindications cannot be accepted by the Julius Zorn GmbH Group.

Side effects

There are no known side effects when using this product correctly. However, if you notice any adverse changes (e.g. irritation of the skin) during the prescribed use, please consult your doctor or medical retailer immediately. If you know that you are allergic to one or more of the product ingredients, please consult your doctor before using the product. The manufacturer does not accept any liability for damage / injuries resulting from improper use or misuse.

In the case of complaints related to this product, such as damage to the knitted fabric or problems with the fit, please contact the medical retailer directly. Only serious incidents that may lead to a major deterioration of the patient's medical condition or to death should be reported to the manufacturer and the relevant authority in the Member State. Serious incidents are defined in Article 2 (65) Regulation (EU) 2017/745 (MDR).

Disposal

You can dispose of your compression garments in the usual way. There are no specific disposal instructions.

Additional information for products containing silver

As silver can deduct heat and electricity, please avoid the direct contact with electric sources. Juzo ScarComfort Silver should not come in contact with metals (e. g. metal clothes hangers) when the textile is wet. The worked in silver could be discoloured thereby. The tanning agents used in the manufacture of the leather / man-made shoe lining may cause discoloration of the stocking and the silver thread in the foot area.

Batik Collection

Juzo compression sleeves in tiedye design are unique products. The colour pattern of the handmade individual items is unique and cannot be reproduced. There may be minor colour changes after the first wash. The characteristic "wrinkled look" originates from the production technique. It is only visible on the unworn product. Due to the special dyeing technique the compression sleeves are slightly longer in unworn condition than products without tie-dye design. However, this does not have any influence on the fitting accuracy and the medical efficacy when worn.

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