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## Purpose

The consistent application of uniform surface pressure by the compression garments has a beneficial effect on scar maturation and development. Therapeutic application of compression garments results in decreased blood circulation and increased tissue elasticity. The aim is to reduce the risk of the formation of hypertrophic scars and/or keloids or achieve a reduction in their appearance. The objective of compression therapy is to flatten the scar, increase tissue elasticity and reduce scar rigidity.

## Notes for application

For compression stockings/pantyhose with open toes, use the Juzo Slippie Donning aid. Your specialist retailer will gladly show you the correct way to apply the product. Wear your compression garments for 23 hours non-stop (or as advised by your doctor). Compression garments are perfect for whole body treatments and they offer a certified UPF 80 UV protection factor in accordance with UV Standard 801. The UV protection only applies to those body regions that are covered with compression garments. Rough, unfiled finger nails or sharp-edged rings etc. may damage the very fine knitted fabric. Therefore please use Juzo Special Gloves. Toe nails and calluses can cause damage, as can defective inner lining of shoes. For stockings with a silicone border, please note: legs should be depilated and free from oils, powder, creams and ointments.

## Compression stockings

1. Turn your compression stocking inside out up to the heel. To do this, pull the stocking over your arm, hold on to the foot section from the inside and pull it out over the heel (Fig. 1).
2. Now stretch the opening with both hands and slip your foot into the stocking. Pull the foot section up to the middle of your foot. It is important that the heel fabric sits right on the heel (Fig. 2).
3. Now, bit by bit, go up the leg from the foot (Fig. 3 and 4).
4. Please do not pull on the knitted fabric or the edges. At the end, the stocking should sit smoothly against your leg without any creases. If you have pulled the stocking up too high, push it back down to the ankle and start again using a little less strength. Tip: before donning, roll over the silicone border (if applicable).

## Compression pantyhose

5. Initially, proceed as in steps 1 – 4 for compression stockings and repeat these steps for the other leg. Pull the pantyhose up to your knees on both legs (Fig. 5 and 6).
6. Now roll the tights up both legs, bit by bit, ensuring the correct fit without creases, particularly at the back of the knee (Fig. 7).
7. At the end, pull the body part up to your waist. The gusset should be positioned vertically in the crotch. Ensure a fit without creases (Fig. 8).

Caution: do not excessively stretch your compression stockings/pantyhose! The compression stockings are moulded to the leg from the bottom to the top in stages.

## Please note:

Only wear compression garments if they have been prescribed by a doctor. Compression garments should only be dispensed by properly trained staff at a medical retailer. Progressing the success of treatment or a change in condition means that the compression garments will have to be remeasured, and, if necessary, new items will need to be made. Your compression garments are largely resistant to oils, ointments, skin moisture and other environmental influences. The durability of the knitted fabric may be affected by certain combinations of the above. This is why it is advisable to have your compression garments checked regularly at a medical retailer. Please contact your specialist retailer in the event of damage. Do not repair the product yourself - this may affect its quality and medical effectiveness. For hygiene reasons, this product is intended for one person only.

## Washing and care instructions

Please note the care instructions on the textile label sewn into your compression garment. Do up any fasteners (if there are any) before washing! We recommend turning compression stockings with silicone borders inside out before washing. Wash compression garments separately when washing for the first time (they may lose colour). Your compression garments should be washed every day on a gentle or delicate wash cycle (40 °C). Tip: a laundry net provides additional protection for the knitted fabric. Please do not use fabric softener. When hand washing, rinse well and do not wring. We recommend the gentle Juzo Special Detergent. You can reduce drying time by placing the compression garment on a thick towel after washing, tightly rolling the towel up and firmly pressing out excess moisture. Then hang up the compression garment loosely. Do not leave it in the towel and do not dry on a radiator or in the sun. Compression garments are not suitable for dry cleaning or tumble drying.

## Material composition

For precise details, please see the textile label that is sewn into your compression garment.

## Storage information and shelf life

Store in a dry place and protect from sunlight. Standard products have a maximum shelf life of 36 months. The compression garment then has a usage period of 6 months. The usage period depends on correct handling (e.g. care, donning and doffing) and normal use of the products. The usage period of the medical product is printed on the box with an hourglass symbol. Custom-made products are for immediate use. Body measurements may change based on the patient's condition and individual circumstances. Custom-made products are therefore designed for a usage period of 6 months. After the 6-month usage period, body measurements must be re-checked by a medical specialist. In the event of significant changes in measurements of the body part that is being treated (e.g. based on therapeutic success, weight change), re-measurement and supply of a new product may be required earlier, before the 6-month usage period comes to an end.

## Indications

Damage to the skin and healing for skin with hypertrophic scars and keloids.

- after injuries / burns / operations
- to improve the scar / wound healing through to completely healing the scar
- to stabilise the surgical outcome

## Absolute contraindications

Peripheral artery disease (PAD) stage III/IV, decompensated heart failure, septic phlebitis, phlegmasia cerulea dolens, acute erysipelas, acute bacterial, viral or allergic inflammation, untreated or advanced skin diseases/contact dermatitis, untreated open wounds, swelling of the extremity of unknown cause, untreated malignant tumours.

## Relative contraindications

Pronounced weeping dermatoses, malignant lymphoedema, peripheral artery disease (PAD) stage I/II, peripheral neuropathy in association with diabetes mellitus, sensory disorders of the extremities, compensated heart failure, skin infections, skin diseases (e.g. blistering dermatoses), intolerance to pressure or the ingredients of the product, rheumatoid arthritis. If there are open wounds in the area to be treated, the affected area must be covered with an appropriate wound dressing or bandage before the compression garment can be applied.

For chronic lymphoedema or lipoedema, we recommend our lymphological products.

Liability for failure to comply with these contraindications cannot be accepted by the Julius Zorn GmbH Group.

## Side effects

There are no known side effects when using this product correctly. However, if you notice any adverse changes (e.g. irritation of the skin) during the prescribed use, please consult your doctor or medical retailer immediately. If you know that you are allergic to one or more of the product ingredients, please consult your doctor before using the product. The manufacturer does not accept any liability for damage / injuries resulting from improper use or misuse.

In the case of complaints related to this product, such as damage to the knitted fabric or problems with the fit, please contact the medical retailer directly. Only serious incidents that may lead to a major deterioration of the patient's medical condition or to death should be reported to the manufacturer and the relevant authority in the Member State. Serious incidents are defined in Article 2 (65) Regulation (EU) 2017/745 (MDR).

## Disposal

You can dispose of your compression garments in the usual way. There are no specific disposal instructions.

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**Julius Zorn GmbH** · Juliusplatz 1 · 86551 Aichach · Germany